



Determinación	Técnica	Matriz	Analitos
Aditivos	HPLC	Aceites transformador/Corcho	Ingramet 39, Nypass SP, Irganox 1076
Aromas	GC	Producto limpieza	2-Methylisoborneol, 3-Methyl-1-butanol, 6-methyl-5-hepten-2-ona, Ethyl lactate, Eucaliptol, Geosmin, Isovalerato de etilo, Limonene, Mentol, Octen-3-ol, Pinene
Bisphenol	GC/MS/MS	Plásticos / Resinas	Bisphenol A, o,p-Bisphenol
Compuestos azoicos	HPLC	Textil/cuero	4-Aminodiphenyl, Benzidine, 4-Chlor-o-toluidine, 2-Naphtylamin, o-Aminazotoluene, 5-Nitro-o-toluidine, p-chloroaniline, 2,4-Diaminoanisoole, 4,4'-Diaminophenylmethane, 3,3'-Dichlorobenzidine, 3,3'-Dimethoxybenzidine, 3,3'-Dimethylbenzidine, 4,4'-Methylenedi-o-toluidine, p-Cresidine, 4,4'-Methylene-bis(2-chloroaniline), 4,4'-Oxydianiline, 4,4'-Thiodianiline, o-Toluidine, 2,4-Toluenediamine, 2,4,5-Trimethylaniline, o-Anisidine, 4-Aminoazobenzene
Filtros Solares	HPLC	Cremas solares	Bis-Ethylhexylphenol methoxyphenyl triazine (TINOSORB), Butylmethoxidibenzoylmethane (PARSOL), Ethyl Hexyl Salicylate, Ethylhexyl Methoxycinnamate (PARSOL MCX), Ethylhexyl triazone, Isopentyl-4-methoxycinnamate, Methylbenzylidene Camphor, Octocrylene, Phenylbenzimidazole Sulfonic Acid (Eusolex 232)
Ftalatos	GC/MS/MS	Textil/cuero/plásticos	Dimethyl Phthalate (DMP), Diethyl Phthalate (DEP), Di-isopropyl Phthalate (DPP), Di-isobutyl Phthalate (DIBP), Dibutyl Phthalate (DBP), Benzyl Butyl Phthalate (BBP), Bis (2-ethylhexy) Adipate (DEHA) + Dioctyl Adipate (DOA), Bis (2-Ethylhexyl) Phthalate (DEHP), Di-n-octyl Phthalate (DNOP), Diisononyl Phthalate (DINP), Di-i-Decyl Phthalate (DIDeP), Diisododecyl Phthalate (DIDP)
Glicoles	GC	Productos de limpieza	Butylene glycol , Diacetone Alcohol, Diethylene glycol monobutyl ether, Diethyleneglycol , Dimethyl ether, Dipropylene glycol methyl ether, Dipropyleneglycol Dibenzoate, Ethylene glycol , Ethylene glycol butyl ether , Glycerol , Neopentyl Glycol, Propyl glycol, Propylene glycol , Propylene glycol butyl ether , Propylene glycol methyl ether
Impurezas	GC	Producto farmacéutico	Impurezas
Isocianato	GC	Adhesivos/selladoras	2,4-TDI, 4,4-MDI



Meclozina	HPLC	Swabs	Meclozina
Nicotina	HPLC	Swabs/Cigarrillos	Nicotina
Parabenos	HPLC	Productos cosmética	Methyl Paraben, Ethyl Paraben, Propyl Paraben
Parafinas	GC	Swabs	Parafinas
PCB's	GC/MS/MS	Aceite de transformador	PCB's
Volátiles residuales	GC	Producto industrial	1,4-dioxane, 1-propanol, 2-propanol, Acetone, Acetonitrile, Benzene, Chloroform, Cyclohexane, Dichloromethane, Diethyl ether, Ethanol, Ethyl acetate, Heptane, Methanol, n-Hexane, Pyridine, Tetrahydrofuran, Toluene, Trichloroethylene